

## Human CD55/DAF Recombinant Protein

Catalog Number: PROTP08174

### Overview

Product Name	Human CD55/DAF Recombinant Protein
Description	Human CD55/DAF Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 35-353aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by BCA assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P08174, 35-353aa; ADPDCGLPPD VPNAQPALEG RTSFPEDTVI TYKCEESFVK IPGEKDSVIC LKGSQWSDIE EFCNRSCEVP TRLNSASLKQ PYITQNYFPV GTVVEYECRP GYRREPLSP KLTLQNLKW STAVEFCKKK SCPNPGEIRN GQIDVPGGIL FGATISFSCN TGYKLFGSTS SFCLISGSSV QWSDPLPECR EIYCPAPPQI DNGIIQGERD HYG YRQSVTY ACNKGFTMIG EHSIYCTVNN DEGEWSGPPP ECRGKSLTSK VPPTVQKPTT VNVPTTEVSP TSQKTTTKTT TPNAQATRST PVSRTTKHFH ETTPNKGSGT TSHHHHHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CD55/DAF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.