

# Der-P1 Der P1 Protein Recombinant Protein

Catalog Number: PROTP08176

## Overview

Product Name	Der-P1 Der P1 Protein Recombinant Protein
Description	The E. coli derived recombinant protein contains the Dermatophagoides pteronyssinus Dust Mite Der P1 protein (a.a. 20-320) and fused to a 6 His Tag at C-terminus, having a total Mw of 34.5kDa, pI 5.6.
Size	100ug, 500ug, 1mg
Tag	
Form	
Source	Escherichia coli
Formulation	60mM NaCl and 50mM Tris-HCl pH 8.0.

## Concentration

## Storage

Der-P1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

## Purity

Protein is greater than 95% pure as determined by 10% SDS-PAGE (coomassie staining).

## Amino Acid Sequence

MSIKTFEEYKKA FNKSYATFEDEEAARKNFLESVKYVQSNGGAINHLSDSLDEFK  
NRFLMSAEAFEHLKTQFDLNAETNACSINGNAPAEIDLRQMRTVTPIRMQGGCGSAWAFS  
GVAATESAYLAYRNQSLDLAEQELVDCASQHGCHGDTIPRGIEYIQHNGVVQESYYRYVA  
REQSCRPNARFGISNYCQIYPPNVNKIREALAQTHSAIAVIIIGIKDLDAFRHYDGRTI  
IQRDNGYQPNYHAVNIVGYSNAQGV DYWIVRNSWDTNWGDNGYGYFAANIDLMMIEEYPYVVILHHHHHH

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Der-P1 Der P1 Protein Recombinant Protein