

Human GUSB Recombinant Protein

Catalog Number: PROTP08236

Overview

Product Name	Human GUSB Recombinant Protein
Description	Human GUSB Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 23-651aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P08236, 23-651aa; LQGGMLYPQE SPSRECKELD GLWSFRADFS DNRRRGFEEQ WYRRPLWESG PTVDMPPVSS FNDISQDWRL RHFVGVVWYE REVILPERWT QDLRTRVVL RIGSAHSYAI VVWNGVDTLEH EGGYLPFEAD ISNLVQVGPL PSRLRITIAI NNTLTPTTLP PGTIQYLTDT SKYPKGYFVQ NTYFDFNYA GLQRSVLLYT TPTTYIDDIT VTTSVEQDSG LVNYQISVKG SNLFKLEVRL LDAENKVVAN GTGTQGQLKV PGVSLWWPYL MHERPAYLYS LEVQLTAQTS LGPVSDFYTL PVGIRTVAVT KSQLFLNGKP FYFHGVNKHE DADIRGKGF D WPLLVKDFNL LRWLGANAFR TSHYPYAEV MQMCDRYGIV VIDECPGVGL ALPQFFNNVS LHHHMQMVEE VVRRDKNHPA VVMWSVANEP ASHLESAGYY LKMVIAHTKS LDPSRPVTFV SNSNYAADKG APYVDVICLN SYSWYHDYH HLELIQLQLA TQFENWYKKY QKPIIQSEYG AETIAGFHQD PPLMFTEEQ KSLLEQYHLG LDQKRRKYV V GELIWNFADF MTEQSPTRVL GNKKGIFTRQ RQPKSAAFL RERYWKIANE TRYPHSVAKS QCLENSLFT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GUSB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.