

Mouse Mdh2 Recombinant Protein

Catalog Number: PROTP08249

Overview

Product Name	Mouse Mdh2 Recombinant Protein
Description	Mouse Mdh2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 25-338aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P08249, 25-338aa; MGSSHHHHHH SGLVPRGSH MAKVAVLGAS GGIGQPLSLL LKNSPLVSRLL TLYDIAHTPG VAADLSHIET
RANVKGYLGP EQLPDCLKGC DVVVIPAGVP RPKGMTRDDL FNTNATIVAT LTAACAQHCP EAMVCIANP VNSTIPITAE VFKKHGVYNP
NKIFGVTTLD IVRANTFVAE LKGLDPARVN VPVIGGHAGK TIPLISQCT PKVDFPQDQL ATLTGRIQEA GTEVVKAKAG AGSATLSMAY
AGARFVFLSV DAMNGKEGVV ECSFVQSKET ECTYFSTPLL LGKKGLEKNL GIGKITPFEE KMIAEAIPEL KASIKKGEDF VKNMK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Mdh2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.