

e.coli DnaJ Recombinant Protein

Catalog Number: PROTP08622

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli DnaJ Recombinant Protein
Description	E. coli DnaJ Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-376aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	25mM Tris-HCl buffer (pH 8.8) containing 5mM DTT, 100mM NaCl, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P08622, 1-376aa; MAKQDYIEIL GVSKTAEHE IRKAYKRLAM KYHPDRNQGD KEAEAKFKEI KEAYEVLTDQ KRAAYDQYG HAAFEQGGMG GGGFGGGADF SDIFGDVFGD IFGGGRGRQR AARGADLRYN MELTLEEAVR GVTKEIRIPT LEECDVCHGS GAKPGTQPQT CPTCHGSGQV QMRQGFFAVQ QTCPHCQGRG TLIKDPCNKC HGHGRVERSK TSLVKIPAGV DTGDRIIRLAG EGEAGEHGAP AGDLYVQVQV KQHPIFEREG NNLYCEVPIN FAMAALGGEI EVPTLDGRVK LKVPGETQTG KLFRMRGKGV KSVRGGAQGD LLCRVVETP VGLNERQKQL LQELQESFGG PTGEHNSPRS KSFFDGVKKF FDDLTR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli DnaJ Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.