

Human KRT16 Recombinant Protein

Catalog Number: PROTP08779

Overview

Product Name	Human KRT16 Recombinant Protein
Description	Human KRT16 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-473aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P08779, 1-473aa; MGSSHHHHHH SGLVPRGSH MGSMTTC SRQ FTSSSSMKGS CGIGGGIGGG SSRISVLAG GSCRAPSTYG
GGLSVSSRFS SGGACGLGGG YGGGFSSSSS FGSGFGGGYG GGLGAGFGGG LGAGFGGGFA GGDGLLVGSE KVTMQNLNDR LASYLDK VRA
LEEANADLEV KIRDWYQRQR PSEIKDYSPY FKTIEDLRNK IIAATIEN AQ PILQIDNARL AADDFRTKYE HELALRQTVE ADVNGLRRVL
DELTARTDL EMQIEGLKEE LAYLRKNHEE EMLALRGQTG GDVNVEMDAA PGVDLSRILN EMRDQYEQMA EKNRRDAETW FLSKTEELNK
EVASNSLVQ SSRSEVELR RVLQGLEIEL QSQLSMKASL ENSLEETKGR YCMQLSIIQG LIGSVEEQLA QLRCEMEQQS QEYQILLDVK
TRLEQEIATY RRLLEGEDAH LSSQQASGQS YSSREVF TSS SSSSRQTRP ILKEQSSSF SQGQSS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human KRT16 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.