

Human recombinant FGF basic/FGF-2 protein

Catalog Number: PROTP09038-8

Overview

Product Name	Human recombinant FGF basic/FGF-2 protein
Description	Human recombinant FGF basic/FGF-2 protein expressed in E.coli with Non-Tagged. Sequence domain: 143-288aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E.coli
Formulation	50mM Tris-HCl buffer (pH 8.0)

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P09038, 143-288aa; MPALPEDGGS GAFPPGHFKD PKRLYCKNGG FFLRIHPDGR VDGVREKSDP HIKLQLQAEE RGVVSIGVVC ANRYLAMKED GRLLASKCVT DECFFFERLE SNNYNTYRSR KYTSWYVALK RTGQYKLGSK TGPQKAILF LPMSAKS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant FGF basic/FGF-2 protein

For Research Use Only. Not for use in diagnostic procedures.