

Mouse P4HB Recombinant Protein

Catalog Number: PROTP09103

Overview

Product Name	Mouse P4HB Recombinant Protein
Description	Mouse P4HB Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 20-509aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P09103, 20-509aa; DALEEDNVL VLKKSNFEEA LAAHKYLLVE FYAPWCGHCK ALAPEYAKAA AKLKAEGSEI RLAKVDATEE
SDLAQQYGVR GYPTIKFFKN GDTASPKEYT AGREADDIVN WLKVRTGPAAT TSLSDTAAAE SLVDSSEVTV IGFFKDVESD SAKQFLLAAE
AIDDIPFGIT SNSGVFSKYQ LDKDGVVLFK KFDEGRNNE GEITKEKLLD FIKHNQLPLV IEFTEQTAPK IFGGEIKTHI LLFLPKSVSD YDGKLSFKR
AAEGFKGKIL FIFIDSHTD NQRILEFFGL KKEECPAVRL ITLEEEMTKY KPESDELTAE KITEFCHRFL EGKIKPHLMS QVPEDWDKQ
PVKVLVGANF EEVAFDEKKN VFVEFYAPWC GHCKQLAPIW DKLGETYKDH ENIIIAKMDS TANEVEAVKV HSFPTLKFFP ASADRTVIDY
NGERTLDGFK KFLESGGQDG AGDDEDLDLE EALEPDM EED DDQKAVKDEL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse P4HB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.