

Human ACAA Recombinant Protein

Catalog Number: PROTP09110

Overview

Product Name	Human ACAA Recombinant Protein
Description	Human ACAA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 27-424aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 1mM DTT, 0.1M NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P09110, 27-424aa; MGSSHHHHHH SGLVPRGSH MLGAPQASA ADVVVVHGRR TAICRAGRGG FKDTTPDELL SAVMTAVLKD
VNLRPEQLGD ICVGNVLQPG AGAIMARIAQ FLSDIPETVP LSTVNRQCSS GLQAVASIAG GIRNGSYDIG MACGVESMSL ADRGNPGNIT
SRLMEKEKAR DCLIPMGITS ENVAERFGIS REKQDTFALA SQQKAARAQS KGCFQAEIVP VTTTVHDDKG TKRSITVTQD EGIRPSTTME
GLAKLKPAFK KDGSTTAGNS SQVSDGAAAI LLARRSKAEE LGLPILGVLR SYAVVGVPD IMGIGPAYAI PVALQKAGLT VSDVDIFEIN
EAFASQAAYC VEKLRLLPPEK VNPLGGAVAL GHPLGCTGAR QVITLLNELK RRGKRAYGVV SMCIGTGMGA AAVFEYPGN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ACAA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.