

Human Glutathione S-Transferase pi 1/GSTP1 Recombinant Protein

Catalog Number: PROTP09211

Overview

Product Name	Human Glutathione S-Transferase pi 1/GSTP1 Recombinant Protein
Description	Human Glutathione S-Transferase pi 1/GSTP1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-210aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.0) containing 30% glycerol, 1mM EDTA, 0.1mM PMSF.

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P09211, 1-210aa; MPPY TVVYFPVRGR CAALRMLLAD QGQSWKEEVV TVETWQEGSL KASCLYGQLP KFQDGLTLY QSNTILRHLG
RTLGLYGKDQ QEAALVDMVN DGVEDLRCKY ISLIYTYEA GKDDYVKALP GQLKPFETLL SQNQGGKTFI VGDQISFADY NLLDLLLLIHE
VLAGCLDAF PLLSAYVGR L SARPKLKAF L ASPEYVNLPI NGNGKQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Glutathione S-Transferase pi 1/GSTP1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.