

Human MMP-10 Recombinant Protein

Catalog Number: PROTP09238

Overview

Product Name	Human MMP-10 Recombinant Protein
Description	Human MMP-10 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 99-476aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P09238, 99-476aa; MGSSHHHHHH SGLVPRGSH MGSFSSFGPM PKWRKTHLTY RIVNYTPDLP RDAVDSAIEK ALKWEEVTP LTF SRLYEGE ADIMISFAVK EHGDFYSFDG PGHSLAHAYP PGPGLYGDIIH FDDDEKWTEG ASGTNLF LVA AHELGHS LGL FHSANTEALM YPLNSFTEL AQFRLSQDDV NGIQSLYGPP PASTEEPLVP TKSVP SGSEM PAKCDPALS F DAISTLRGEY LFFKDRYFWR RSHWNPEPEF HLISAFWPSL PSYLDAAYEV NSRDTVFIFK GNEFWAIRGN EVQAGYPRGI HTLGF PPTIR KIDAAVSDKE KKKTYFFAAD KYWRFDENSQ SMEQGFPRLI ADDFP GVEPK VDAVLQAFGF FYFFSGSSQF EFDPNARMVT HILKSNSWLH C

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MMP-10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.