

Mouse Alkaline Phosphatase/ALPL Recombinant Protein

Catalog Number: PROTP09242

Overview

Product Name	Mouse Alkaline Phosphatase/ALPL Recombinant Protein
Description	Mouse Alkaline Phosphatase/ALPL Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 19-503aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P09242, 19-503aa; VPEKERDPSY WRQQAQETLK NALKLQKLNT NVAKNVIMFL GDGMGVSTVT AARILKGQLH HNTGEETRL E
MDKFPFVALS KTYNTNAQVP DSAGTATAYL CGVKANEGTV GVSAATERTR CNTTQGNEVT SILRWAKDAG KSVGIVTTTR VNHATPSAAY
AHSADRDWYS DNEMPPEALS QGCKDIAYQL MHNKIDIDVI MGGGRKYMYP KNRTDVEYEL DEKARGTRLD GLDLISIWKS FKPRHKHSHY
VWNRTELLAL DPSRVDYLLG LFEPGDMQYE LNRNLTDP S LSEMVEVALR ILTKNLKGF LLVEGGRIDH GHHEGKAKQA LHEAVEMDQA
IGKAGAMTSQ KDTLTVVTAD HSHVFTFGGY TPRGNSIFGL APMVSDTDKK PFTAILYGNG PGYKVVDGER ENVSMVDYAH NNYQAQSAVP
LRHETHGGED VAVFAKGPMA HLLHGVHEQN YIPVHMAYAS CIGANLDHCA WAGSG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Alkaline Phosphatase/ALPL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.