

Mouse Pgk1 Recombinant Protein

Catalog Number: PROTP09411

Overview

Product Name	Mouse Pgk1 Recombinant Protein
Description	Mouse Pgk1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-417aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P09411, 1-417aa; MGSSHHHHHH SGLVPRGSH MGSHMSLSNK LTLDKLDVKG KRVVMRVDFN VPMKNNQITN NQRKAAVPS IKFCLDNGAK SVVLMShLGR PDGVPMPDKY SLEPVAAELK SLLGKDVLF LKDCVGPVEVEN ACANPAAGTV ILLENLRFHV EEEGKGDAS GNKVKAEPAK IDAFRASLSK LGDVYVND AF GTAHRAHSSM VGVNLPQKAG GFLMKKELNY FAKALES PER PFLAILGGAK VADKIQLINN MLDKVNEMII GGGMAFTFLK VLNNMEIGTS LYDEEGAKIV KDLMSKAEKN GVKITLPVDF VTADKFDENA KTGQATVASG IPAGWMGLDC GTESSKKYAE AVGRAKQIVW NGPVGVFEWE AFARGTKSLM DEVVKATSRG CITIIGGGDT ATCCAkwNTE DKVSHVSTGG GASLELLE GK VLPGVDA LSN V

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Pgk1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.