

## Human Annexin A4 Recombinant Protein

Catalog Number: PROTP09525

### Overview

Product Name	Human Annexin A4 Recombinant Protein
Description	Human Annexin A4 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-321aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol 0.2M NaCl

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P09525, 1-321aa; MAMATKGGTV KAASGFNAME DAQTLRKAMK GLGTDEDAIL SVLAYRNTAQ RQEIRTAYKS TIGRDLIDDL  
KSELSGNFEQ VIVGMMTPTV LYDVQELRRA MKGAGTDEGC LIEILASRTP EEIRRISQTY QQYGRSLED DIRSDTSFMF QRVLVLSL SAG  
GRDEGNYLDD ALVRQDAQDL YEAGEKKWGT DEVKFLTVLC SRNRNHLLHV FDEYKRISQK DIEQSIKSET SGSFEDALLA IVKCMRNKSA  
YFAEKLYKSM KGLGTDDNTL IRVMVSRAEI DMLDIRAHFK RLYGKSLYSF IKGDTSGDYR KVLLVLCGGD D

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Annexin A4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.