

Human PARP Recombinant Protein

Catalog Number: PROTP09874

Overview

Product Name	Human PARP Recombinant Protein
Description	Human PARP Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 662-1014aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P09874, 662-1014aa; MKSKLPKPVQ DLIKMIFDVE SMKKAMVEYE IDLQKMPLGK LSKRQIQAA Y SILSEVQQAV SQGSSDSQIL
DLSNRFYTLI PHDFGMKKPP LLNADSVQA KAEMLDNLLD IEVAYSLLRG GSDDSSKDPI DVNYEKLKTD IKVVDRDSEE AEIIRKYVKN
THATTHNAYD LEVIDIFKIE REGECQRYKP FKQLHNRLL WHGSRTTNFA GILSGLRIA PPEAPVTGYM FGKGIYFADM VSKSANYCHT
SQGDPIGLIL LGEVALGNMY ELKHASHISK LPKKGKHSVKG LGKTTDPDSA NISLDGVDVP LGTGISSGVN DTSLLYNEYI VYDIAQVNLK YLLKLFNFK
TSLW

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PARP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.