

## Human UCH-L1/PGP9.5 Recombinant Protein

Catalog Number: PROTP09936

### Overview

Product Name	Human UCH-L1/PGP9.5 Recombinant Protein
Description	Human UCH-L1/PGP9.5 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-223aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM EDTA

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P09936, 1-223aa; MQLKPMEINP EMLNKVLSRL GVAGQWRFVD VLGLEESLG SVPAPACALL LLFPLTAQHE NFRKKQIEEL  
KGQEVSPKVY FMKQTIGNSC GTIGLIHAVA NNQDKLGFED GSVLKQFLSE TEKMSPEDRA KCFEKNEAIQ AAHDAVAQEG QCRVDDKVN  
HFILFNNVDG HLYELDGRMP FPNVNHGASSE DTLLKDAAKV CREFTEREQG EVRFSVAVALC KAA

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human UCH-L1/PGP9.5 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.