

e.coli ackA Recombinant Protein

Catalog Number: PROTP0A6A3

Overview

Product Name	e.coli ackA Recombinant Protein
Description	E. coli ackA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-400aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A6A3, 1-400aa; MGSSHHHHHH SGLVPRGSH MGSMSSKLVL VLNCGSSSLK FAIDAVNGE EYLSGLAECF HLPEARIKWK
MDGNKQEAAL GAGAAHSEAL NFIVNTILAQ KPESLAQLTA IGHRIVHGGE KYTSSVIDE SVIQGIKDA SFAPLHNPALH LIGIEEALKS
FPQLKDKNVA VFDTAFHQTM PEESYLALP YNLYKEHGIR RYGAHGTSHF YVTQEAAMKML NKPVEELNII TCHLGNGGSV SAIRNGKCVD
TSMGLTPLEG LVMGTRSGDI DPAIIFHLHD TLGMSVDAIN KLLTKESGLL GLTEVTSDCR YVEDNYATKE DAKRAMDVYC HRLAKYIGAY
TALMDGRLDA VVFTGGIGEN AAMVRELSLG KLGVLGFEVD HERNLAARFG KSGFINKEGT RPAVV IPTNE ELVIAQDASR LTA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli ackA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.