

e.coli glpK Recombinant Protein

Catalog Number: PROTP0A6F3

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli glpK Recombinant Protein
Description	E. coli glpK Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-502aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A6F3, 1-502aa; MGSSHHHHHH SGLVPRGSH MGSMTTEKKYI VALDQGTSS RAVVMDHDAN IISVSQREFE QIYPKPGWVE
HDPMEIWATQ SSTLVEVLAK ADISSDQIAA IGITNQRETT IVWEKETGKP IYNAIVWQCR RTAEICEHLK RDGLEDYIRS NTGLVIDPYF
SGTKVKWILD HVEGSRERAR RGELLFGTVD TWLIWKMTQG RVHVTDYTNA SRTMLFNIHT LDWDDKMLEV LDIPREMLPE VRRSSEVYGQ
TNIGGKGGTR IPISGIAGDQ QAALFGQLCV KEGMAKN TYG TGCFMLMNTG EKAVKSENG LTTIACGPTG EVNYALEGAV FMAGASIQWL
RDEMKLINDA YDSEYFATKV QNTNGVYVVP AFTGLGAPYW DPYARGAIFG LTRGVNANH IIRATLESIA YQTRDVL EAMQ ADSGIRLHAL

RVDGGAVANN FLMQFQSDIL GTRVERPEVR EVTALGAAYL AGLAVGFWQN LDELQEKAVI EREFRPGIET TERNYRYAGW KKAVKRAMAW
EEHDE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli glpK Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.