

e.coli coaA Recombinant Protein

Catalog Number: PROTP0A6I3

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli coaA Recombinant Protein
Description	E. coli coaA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-316aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 10% glycerol, 200mM NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A6I3, 1-316aa; MGSSHHHHHH SSGLVPRGSH MGS HMSIQE TLMPYLQFD RNQWAALRDS VPMTLSEDEI ARLKGINEDL SLEEVAEIYL
PLSRLLNFYI SSNLRRQAVL EQFLGTNGQR IPYIISIAGS VAVGKSTTAR VLQALLSRWP EHRRVELITT DGFLHPNQVL KERGLMKKKG
FPSYDMHRL VKFVSDLKSG VPNTAPVYS HLIYDVIPDG DKTVVQPDIL ILEGLNVLQS GMDYPHDPHH VFVSDFVDFS IYVDAPEDLL
QTWYINRFLK FREGAFTDPD SYFHNYAKLT KEEAIKTAMT LWKEINWNLN KQNILPTRER ASLILTKSAN HAVEEVRLRK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli coaA Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.