

e.coli nanA Recombinant Protein

Catalog Number: PROTP0A6L4

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli nanA Recombinant Protein
Description	E. coli nanA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-297aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A6L4, 1-297aa; MGSSHHHHHH SGLVPRGSH MATNLRGVMA ALLTPFDQQQ ALDKASLRRL VQFNIQQGID GLYVGGSTGE
AFVQSLSERE QVLEIVAEAA KGKIKLIAHV GCVSTAESQQ LAASAKRYGF DAVSAVTPFY YPFSFEEHCD HYRAIIDSAD GLPMVVYNIP
ALSGVKLTLD QINTLVTLPG VGALKQTSKD LYQMEQIRRE HPDLVLYNGY DEIFASGLLA GADGGIGSTY NIMGWRYQGI VKALKEGDIQ
TAQKLQTECN KVIDLLIKTG VFRGLKTVLH YMDVVSPLC RKPFGPVDEK YLPELKALAQ QLMQERG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli nanA Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.