

## e.coli Dnak Recombinant Protein

Catalog Number: PROTP0A6Y8

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                           |
|--------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name | e.coli Dnak Recombinant Protein                                                                                           |
| Description  | E. coli Dnak Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 385-546aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug                                                                                                        |
| Tag          | Non-Tagged                                                                                                                |
| Form         | Liquid                                                                                                                    |
| Source       | E. coli                                                                                                                   |
| Formulation  | 25mM Tris-HCl buffer (pH 7.5) containing 2mM beta-mercaptoethanol, 1mM EDTA                                               |

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

E.coli, P0A6Y8, 385-546aa; MDVKDVLLLD VTPLSLGIET MGGVMTTLIA KNTTIPTKHS QVFSTAEDNQ SAVTIHVLQG ERKRAADNKS  
LGQFNLDGIN PAPRGMPQIE VTFDIDADGI LHVSAKDKNS GKEQKITIKA SSGLNEDEIQ KMRVDAEANA EADRKFEELV QTRNQGDHLL HST

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### [e.coli Dnak Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.