

E.coli recombinant rnpA protein

Catalog Number: PROTP0A7Y8

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	E.coli recombinant rnpA protein
Description	E.coli recombinant rnpA protein expressed in E.coli with Non-Tagged. Sequence domain: 1-119aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E.coli
Formulation	Liquid. Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A7Y8, 1-119aa; MVKLAFPREL RLLTPSQFTF VFQQPQRAGT PQITILGRIN SLGHPRIGLT VAKKNVRRAH ERNRIKRLTR ESFRLRQHEL PAMDFVVVAK KGVADLDNRA LSEALEKLWR RHCRLARGS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[E.coli recombinant rnpA protein](#)

For Research Use Only. Not for use in diagnostic procedures.