

e.coli tnaA Recombinant Protein

Catalog Number: PROTP0A853

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli tnaA Recombinant Protein
Description	E. coli tnaA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-471aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A853, 1-471aa; MENFKHL PEPFRIRVIE PVKRTTRAYR EEAIKSGMN PFLDSEDVF IDLLTDSGTG AVTQSMQAAM MRGDEAYSGS RSYALAESV KNIFGYQYTI PTHQGRGAEQ IYIPVLIKKR EQEKGLDRSK MVAFSNYFFD TTQGHSQLING CTVRNVYIKE AFDTGVRDYDF KGNFDLEGLE RGIEEVGPNN VPYIVATITS NSAGGQPVSL ANLKAMYSIA KKYDIPVVM D SARFAENAYF IKQREAEYKD WTIEQITRET YKYADMLAMS AKKDAMVPMG GLLCMKDDSF FDVYTECRTL CVVQEGFPTY GGLEGGAMER LAVGLYDGMN LDWLAYRIAQ VQYLVDGLEE IGVVCQQAGG HAAFVDAGKL LPHIPADQFP AQALACELYK VAGIRAVEIG SFLGRDPKT GKQLPCPAEL LRLTIPRATY TQTHMDFIIE

AFKHVKENAA NIKGLTFTYE PKVLRHFTAK LKEV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli tnaA Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.