

e.coli TRXB Recombinant Protein

Catalog Number: PROTP0A9P4

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli TRXB Recombinant Protein
Description	E. coli TRXB Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-321aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

E.coli, P0A9P4, 1-321aa; MGTTKHSKLL ILGSGPAGYT AAVYAARANL QPVLITGMEK GGQLTTTTEV ENWPGDPNDL TGPLLMERMH
EHATKFETEI IFDHINKVDL QNRPFRNLNGD NGEYTCDALI IATGASARYL GLPSEEAFCG RGVSACATCD GFFYRNQKVA VIGGGNTAVE
EALYLSNIAS EVHLIHRDGD FRAEKILIKR LMDKVENGNI ILHTNRTLEE VTGDQMGVTG VRLRDTQNSD NIESLDVAGL FVAIGHSPNT AIFEGQLELE
NGYIKVQSGI HGNAQTQSIP GVFAAGDVMD HIYRQAITS A GTGCMAALDA ERYLDGLADA K

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli TRXB Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.