

## e.coli G6PD Recombinant Protein

Catalog Number: PROTP0AC53

### Overview

Product Name	e.coli G6PD Recombinant Protein
Description	E. coli G6PD Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-491aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 0.1mM PMSF, 2mM EDTA, 0.5mM DTT

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

E.coli, P0AC53, 1-491aa; MAVTQTAQAC DLVIFGAKGD LARRKLLPSL YQLEKAGQLN PDTRIIGVGR ADWDKAAAYTK VVREALETFM KETIDEGPWD TLSARLDFCN LDVNDTAAFS RLGAMLDQKN RITINYFAMP PSTFGAICKG LGEAKLNAKP ARVVMKPLG TSLATSQEIN DQVGEYFEEC QVYRIDHYLG KETVLNLLAL RFANSLFVNN WDNRTIDHVE ITVAEEVGIE GRWGYFDKAG QMRDMIQNHLLQILCMIAMS PPSDLSADSI RDEKVKVLKS LRRIDRSNVR EKTVRGQYTA GFAQGKKVPG YLEEEGANKS SNTETFVAIR VDIDNWRWAG VPFYLRGKR LPTKCSEVVV YFKTPELNLF KESWQDLPQN KLTIRLQDPE GVDIQVLNKKV PGLDHHKHLQ ITKLDLSYSE TFNQTHLADA YERLLETMR GIQALFVRRD EVEEAWKWVD SITEAWAMDN DAPKPYQAGT WGPVASVAMI TRDGRSWNEF E

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli G6PD Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.