

e.coli pykF Recombinant Protein

Catalog Number: PROTP0AD61

Overview

Product Name	e.coli pykF Recombinant Protein
Description	E. coli pykF Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-470aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0AD61, 1-470aa; MGSSHHHHHH SGLVPRGSH MGSHMKKTKI VCTIGPKTES EEMLAKMLDA GMNVMRLNFS HGDYAEHGQR
IQNLRNVMSK TGKTAAILLD TKGPEIRTMK LEGGNDVSLK AGQTFFTTDD KSVIGNSEMV AVTYEGFTTD LSVGNTVLVD DGLIGMEVTA
IEGNKVICVK LNNGLGENK GVNLPGVSA LPALAEKDKQ DLIFGCEQGV DFVAASFIRK RSDVIEIREH LKAHGGENIH IISKIENQEG LNNFDEILEA
SDGIMVARGD LGVEIPVEEV IFAQKMMIEK CIRARKVVIT ATQMLDSMIK NPRPTRAEG DVANAILDGT DAVMLSGESA KGKYPLEAVS
IMATICERTD RVMNSRLEFN NDNRLRITE AVCRGAVETA EKLDAPLIVV ATQGGKSARA VRKYFPDATI LALTNEKTA HQLVLSKGVV
PQLVKEITST DDFYRLGKEL ALQSGLAHKG DVVVMVSGAL VPSGTTNTAS VHVL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli pykF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.