

## e.coli Methionine Aminopeptidase/METAP Recombinant Protein

Catalog Number: PROTP0AE18

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                                                    |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | e.coli Methionine Aminopeptidase/METAP Recombinant Protein                                                                                         |
| Description  | E. coli Methionine Aminopeptidase/METAP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-264aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug                                                                                                                                 |
| Tag          | His-Tag                                                                                                                                            |
| Form         | Liquid                                                                                                                                             |
| Source       | E. coli                                                                                                                                            |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT                                                                                     |

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

E.coli, P0AE18, 1-264aa; MGSSHHHHHH SSGLVPRGSH MAISIKTPED IEKMRVAGRL AAEVLEMIEP YVKPGVSTGE LDRICNDYIV NEQHAVSACL GYHGYPKSV C ISINEVVCHG IPDDAKLLKD GDIVNIDVTV IKDGFHGDTS KMFIVGKPTI MGERLCRITQ ESLYLALRMV KPGINLREIG AAIQKFVEAE GFSVVREYCG HGIGRGFHEE PQVLHYDSRE TNVVLKPGMT FTIEPMVNAG KKEIRTMKDG WTVKTKDRSL SAQYEHTIVV TDNGCEILTL RKDDTIPAI SHDE

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli Methionine Aminopeptidase/METAP Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.