

e.coli Dsbc Recombinant Protein

Catalog Number: PROTP0AEG6

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli Dsbc Recombinant Protein
Description	E. coli Dsbc Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 21-236aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 2mM EDTA

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0AEG6, 21-236aa; MDDAAIQQTL AKMGIKSSDI QPAPVAGMKT VLTNSGVLYI TDDGKHIIQG PMYDVSGTAP VNVTNKMLLK
QLNALEKEMI VYKAPQEKHV ITVFTDITCG YCHKLHEQMA DYNALGITVR YLAFPRQGLD SDAEKEMKAI WCAKDKNKAF DDVMAGKSVA
PASCDVDIAD HYALGVQLGV SGTPAVVLSN GTLVPGYQPP KEMKEFLDEH QKMTSGK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli Dsbc Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.