

e.coli SecB Recombinant Protein

Catalog Number: PROTP0AG86

Overview

Product Name	e.coli SecB Recombinant Protein
Description	E. coli SecB Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-155aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

1mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P0AG86, 1-155aa; MSEQNNTMT FQIQRIYTKD ISFEAPNAPH VFQKDWQPEV KLDLDTASSQ LADDVYEVVL RVTVTASLGE
ETAFLCEVQQ GGIFSIAGIE GTQMAHCLGA YCPNILFPYA RECITSMVSR GTFPQLNLAP VNFDAFMNY LQQQAGEGTE EHQDA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli SecB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.