

Human recombinant ARSM2 protein

Catalog Number: PROTP0C7Q2-1

Overview

Product Name	Human recombinant ARSM2 protein
Description	Human recombinant ARSM2 protein expressed in E.coli with His-Tag. Sequence domain: 1-107aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E.coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing , 10% glycerol, 0.4M urea

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P0C7Q2, 1-107aa; MGSSHHHHHH SSGLVPRGSH MGSMRLRYPG PMVTEAEGKG GPEMASLSSS VVPVSFISTL RESVLDPGVG
GEGASDKQRS KLSLSHSMIP AAKIHTELCL PAFFSPAGTQ RRFQQPQHHL TLSIIHTAAR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant ARSM2 protein

For Research Use Only. Not for use in diagnostic procedures.