

Mouse Annexin A1 Recombinant Protein

Catalog Number: PROTP10107

Overview

Product Name	Mouse Annexin A1 Recombinant Protein
Description	Mouse Annexin A1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-346aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P10107, 1-346aa; MAMVSEF LKQARFLENQ EQEYVQAVKS YKGGPGSAVS PYPSFNVSSD VAALHKAIMV KGVDEATIID ILTKRTNAQR QQIKAAYLQE NGKPLDEVLR KALTGHLEEV VLAMLKTPAQ FDADELRGAM KGLGTDEDL IEILTTRSNE QIREINRVYR EELKRDLAKD ITSDTSGDFR KALLALAKGD RCQDLSVNQD LADTDARALY EAGERRKGTD VNVFTTILTS RSFPHLRRVF QNYGKYSQHD MNKALDLELK GDIEKCLTTI VKCATSTPAF FAEKLYEAMK GAGTRHKALI RIMVSRSEID MNEIKVFYQK KYGISLCQAI LDETKGDYEK ILVALCGGN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Annexin A1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.