

Mouse CCL2/JE/MCP-1 Recombinant Protein

Catalog Number: PROTP10148

Overview

| | |
|--------------|---|
| Product Name | Mouse CCL2/JE/MCP-1 Recombinant Protein |
| Description | Mouse CCL2/JE/MCP-1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-148aa. Application(s): SDS-PAGE. |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | Phosphate-Buffered Saline (pH 7.4) 10% glycerol |

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Mouse, P10148, 24-148aa: QPDAVNAPL TCCYSFTSKM IPMSRLESYK RITSSRCPKE AVVFVTKLKR EVCADPKKEW VQTYIKNLDR NQMRSEPTTL FKTASALRSS APLNVKLRK SEANASTTFS TTSSTSVGV TSVTVN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse CCL2/JE/MCP-1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.