

Human Tau Recombinant Protein

Catalog Number: PROTP10636

Overview

Product Name	Human Tau Recombinant Protein
Description	Human Tau Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-352aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 1mM DTT, 10% glycerol

Concentration

0.5mg/ml (determined by BCA assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P10636, 1-352aa; MAEPRQEFEV MEDHAGTYGL GDRKDQGGYT MHQDQEGD TD AGLKAE EAGI GDTPSLEDEA AGHVTQARMV
SKSKDGTGSD DKKAKGADGK TKIATPRGAA PPGQKQANA TRIPAKTPPA PKTPPSSGEP PKSGDRSGYS SPGSPGTPGS RSRTPSLPTP
PTREPKKVAV VRTPPKSPSS AKSRLQTAPV PMPDLKNVKS KIGSTENLKH QPGGGKVQIV YKPVDLSKVT SKCGSLGNIH HKPGGGQVEV
KSEKLDLFKDR VQSKIGSLDN ITHVPGGGNK KIETHKLTFR ENAKAKTDHG AEIVYKSPVV SGGTSPRHLS NVSSTGSIDM VDSPQLATLA
DEVASLAKQ GL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Tau Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.