

Human recombinant HAPLN1 protein

Catalog Number: PROTP10915-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant HAPLN1 protein
Description	Human recombinant HAPLN1 protein expressed in HEK293 with His-Tag. Sequence domain: 16-354aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P10915, 16-354aa; DHLSDN YTL DHDRAIHIQAE NGPHLLVEAE QAKVFSHRGG NVTLPCKFYR DPTAFGSGIH KIRIKWTKLT SDYLKEVDVF VSMGYHKKTY GGYQGRVFLK GGS DSDASLV ITDLTLEDYG RYKCEVIEGL EDDTVVVALD LQGVVFPYFP RLGRYNLNFH EAQQACLDQD AVIASFDQLY DAWRGGLDWC NAGWLS DGSV QYPITKPREP CGGQNTVPGV RNYGFWDKDK SRYDVCFTS NFNGRFYYLI HPTKLT YDEA VQACLNDGAQ IAKVGQIFAA WKILGYDRCD AGWLADGSVR YPISRPRRRC SPTEAAVR FV GFPDKKHKLY GVYCFRAYN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant HAPLN1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.