

Human HSPA8/HSC71 Recombinant Protein

Catalog Number: PROTP11142

Overview

Product Name	Human HSPA8/HSC71 Recombinant Protein
Description	Human HSPA8/HSC71 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-646aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P11142, 1-646aa; MSKGPVAVGID LGTTYSCVGV FQHGKVEIIA NDQGNRTTPS YVAFTDTERL IGDAAKNQVA MNPTNTVFDA
KRLIGRRFDD AVVQSDMKHW PFMVVNDAGR PKVQVEYKGE TKSFYPEEVS SMVLTKMKEI AEAYLGKTVT NAVVTVPAYF NDSQRQATKD
AGTIAGLNVL RIINEPTAAA IAYGLDKKVG AERNVLIFDL GGGTFDVSIL TIEDGIFEVK STAGDTHLGG EDFDNRMVNH FIAEFKRKHK
KDISENKRAV RRLRTACERA KRTLSSTQA SIEIDSLYEG IDFYTSITRA RFEELNADLF RGTLDPVEKA LRDAKLDKSQ IHDIVLVGGS TRIPKIQKLL
QDFFNKELN KSINPDEAVA YGAAVQAAIL SGDKSENVQD LLLLDVTPLS LGIETAGGVM TVLIKRNTTI PTKQTQTFTT YSDNQPGVLI
QVYEGERAMT KDNLLGKFE LTGIPPAPRG VPQIEVTFDI DANGILNVA VDKSTGKENK ITITNDKGRL SKEDIERMVQ EAEKYKAEDE
KQRDKVSSKN SLESYAFNMK ATVEDEKLQG KINDEDKQKI LDKCNEIINW LDKNQTAKE EFEHQQKELE KVCNPIITKL YQSAGGMPGG
MPGGFPGGGA PPSGGASSGP TIEEVD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HSPA8/HSC71 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.