

## Mouse recombinant Midkine protein, AF

Catalog Number: PROTP12025-2

### Overview

Product Name	Mouse recombinant Midkine protein, AF
Description	Midkine, also known as neurite growth-promoting factor 2 (NEGF2) is a member of a small family of secreted Growth Factors, furthermore high expression in embryo, and limb E14.5. Mouse midkine shares 87% sequence identity with human midkine. Midkine is a 13.5 kDa protein containing 140 amino acids, which promotes angiogenesis, cell growth, migration, and gene expression of different cell types probably via a multiprotein receptor complex consisting of several molecules.
Size	5ug,20ug,100ug
Tag	His Tag (C-term)
Form	Lyophilized
Source	Escherichia coli
Formulation	The protein was lyophilized from a 0.2 um filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

### Storage

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

### Purity

>98% as determined by SDS-PAGE.

### Molecular weight

The protein has a calculated MW of 14.02 kDa. The protein migrates about 17 kDa under reducing condition (SDS-PAGE analysis).

### Activity

Testing in process

### Endotoxin

<0.1 EU per 1 ug of the protein by the LAL method.

### Amino Acid Sequence

MKKKEKVKKGSECSEWTWGPCTPSSKDCGMGFREGTCGAQTQRVHCKVPCNWKKEFGADCKYKFESWGACDGGSTGTKARQGTLKKARYNAQCQ  
ETIRVTKPCTSKTKSKTKAKKGGKGD with polyhistidine tag at the C-terminus.

### Reconstitution

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H<sub>2</sub>O to a concentration not less than 100 ug/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant Midkine protein, AF

For Research Use Only. Not for use in diagnostic procedures.