

Mouse Pentraxin2/SAP Recombinant Protein

Catalog Number: PROTP12246

Overview

Product Name	Mouse Pentraxin2/SAP Recombinant Protein
Description	Mouse Pentraxin2/SAP Recombinant Protein expressed in HEK293 with His-tag. Sequence domain: 21-224aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P12246, 21-224aa; QTDLKRKVFV FPRESETDHV KLIPHLEKPL QNFTLCFR TYSDLRSQSLF SYSVKGRDNE LLIYKEKVG EYSLYIGQSKV TVRGMEEYLS PVHLCCTWES SSGIVEFWVN GKPWWKLSLQ REYTVKAPPS IVLGQE QDNYGGGFQRSQSF VGEFSDLYMW DYVLT PQDIL FVYRDSPVNP NILNWQALNY EINGYVVIRP RVWD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Pentraxin2/SAP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.