

e.coli ung Recombinant Protein

Catalog Number: PROTP12295

Overview

Product Name	e.coli ung Recombinant Protein
Description	E. coli ung Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-229aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

E.coli, P12295, 1-229aa; MANELTW HDVLAEEKQQ PYFLNLTQTV ASERQSGVTI YPPQKDFVNA FRFTELGDVK VVILGQDPYH GPGQAHGLAF SVRPGIAIPP SLLNMYKELE NTIPGFTRPN HGYLESWARQ GVLLNLTVLT VRAGQAHSHA SLGWETFTDK VISLINQHRE GVVFLWGSN AQKKGAIIDK QRHHVLKAPH PSPLSAHRGF FGCNHFVLN QWLEQRGETP IDWMPVLP AE SE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli ung Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.