

## Human recombinant Fc epsilon RI alpha protein

Catalog Number: PROTP12319-3

### Overview

Product Name	Human recombinant Fc epsilon RI alpha protein
Description	Human recombinant Fc epsilon RI alpha protein expressed in HEK293 with His-Tag. Sequence domain: 26-205aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P12319, 26-205aa; VPQKPKVSLN PPWNRIFKGE NVT LTCNGNN FFEVSSTKWF HNGSLSEETN SSLNIVNAKF EDSGEYKCQH  
QQVNESEPVY LEVFSDWLLL QASAEVVM EG QPLFLRCHGW RNWDVYKVIY YKDGEALKYW YENHNISITN ATVEDSGTY Y CTGKVVQLDY  
ESEPLNITVI KAPREKYWLQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant Fc epsilon RI alpha protein

For Research Use Only. Not for use in diagnostic procedures.