

hiv-1 Gag p24 Recombinant Protein

Catalog Number: PROTP12493

Overview

Product Name	hiv-1 Gag p24 Recombinant Protein
Description	Hiv-1 Gag p24 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 155-321aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 0.1mM PMSF, 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 80% by SDS-PAGE

Amino Acid Sequence

HIV-1, P12493, 155-321aa; WVKVVEEKA FSPEVIMFS ALSEGATPQD LNTMLNTVGG HQAAMQMLKE TINEEAAEWD RLHPVHAGPI
APGQMREPRG SDIAGTTSTL QEQIGWMTHN PPIPVGEIYK RWIILGLNKI VRMYSPTSIL DIRQGPKEPF RDYVDRFYKT LRAEQASQEV
KNWMTETL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



hiv-1 Gag p24 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.