

Human CKMT1A Recombinant Protein

Catalog Number: PROTP12532-1

Overview

Product Name	Human CKMT1A Recombinant Protein
Description	Human CKMT1A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 40-417aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P12532, 40-417aa; MGSSHHHHHH SGLVPRGSH MGSHTMASERR RLYPPSAEYP DLRKHNNCMA SHLTPAVYAR LCDKTTPTGW
TLDQCIQTGV DNPQHPIKT VGMVAGDEET YEVFADLFDV VIQERHNGYD PRTMKHTTDL DASKIRSGYF DERYVLSSRV RTGRSIRGLS
LPPACTRAER REVERVVVDA LSGLKGDLAG RYYRLSEMTE AEQQQLIDDH FLFDKPVSP LTAAGMARDW PDARGIWHNN EKSFLIWWNE
EDHTRVISM KGGNMKRVFE RFCRGLKEVE RLIQERGWEF MWNERLGYIL TCPSNLGTGL RAGVHIKPL LSKDSRFPKI LENLRLQKRG
TGGVDTAATG GVFDISNLDR LGKSEVELVQ LVIDGVNYLI DCERRLERGQ DIRIPTPVIH TKH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CKMT1A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.