

## Human BCKDHA Recombinant Protein

Catalog Number: PROTP12694

### Overview

Product Name	Human BCKDHA Recombinant Protein
Description	Human BCKDHA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 46-445aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 30% glycerol, 0.2M NaCl

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 80% by SDS-PAGE

### Amino Acid Sequence

Human, P12694, 46-445aa; MSSLDDKPQF PGASAEFIDK LEFIQPNVIS GIPIYRVMDR QGQIINPSED PHLPKKVLK LYKSMTLLNT MDRILYESQR  
QGRISFYMTN YGEEGTHVGS AAALDNTDLV FGQYREAGVL MYRDYPLELF MAQCYGNISD LGKGRQMPVH YGCKERHFVT ISSPLATQIP  
QAVGAAAYAAK RANANRVVIC YFGEGAASEG DAHAGFNFAA TLECPIIFFC RNNGYAISTP TSEQYRGDGI AARGPGYGIM SIRVDGNDVF  
AVYNATKEAR RRAVAENQPF LIEAMTYRIG HHSTSDSSA YRSVDEVNYW DKQDHPISRL RHYLLSQGWV DEEQEKAWRK QSRRKVM EAF  
EQAERKPKPN PNLLFSDVYQ EMPAQLRKQQ ESLARHLQTY GEHYPLDHF D K

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human BCKDHA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.