

Mouse FABP1/L-FABP Recombinant Protein

Catalog Number: PROTP12710

Overview

Product Name	Mouse FABP1/L-FABP Recombinant Protein
Description	Mouse FABP1/L-FABP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-127aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P12710, 1-127aa; MGSSHHHHHH SGLVPRGSH MGSMNFGSKY QLQSQENFEP FMKAIGLPED LIQKGKDIKG VSEIVHEGKK IKLTITYGPK VVRNEFTLGE ECELETMTGE KVKAVVKLEG DNKMVTTFKG IKSVELNGD TITNTMTLGD IVYKRYSKRI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse FABP1/L-FABP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.