

e.coli uDP Recombinant Protein

Catalog Number: PROTP12758

Overview

Product Name	e.coli uDP Recombinant Protein
Description	E. coli uDP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-253aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT, 50mM NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P12758, 1-253aa; MGSSHHHHHH SSGLVPRGSH MSKSDVFHLG LTKNDLQGAT LAIVPGDPDR VEKIAALMDK PVKLASHREF
TTWRAELDGK PVIVCSTGIG GPSTSI AVEE LAQLGIRTF L RIGTTGAIQP HINVGDLVLT TASVRLDGAS LHFAPLEFPA VADFECTTAL VEA AKSIGAT
THVGV TASSD TFYPGQERYD TYSGRVVRHF KGSMEEWQAM GVMNYEMESA TLLTMCASQG LRAMVAGVI VNRTQQEIPN AETMKQTESH
AVKIVVEAAR RLL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli uDP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.