

Human recombinant NCAM-1/CD56 protein

Catalog Number: PROTP13591

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant NCAM-1/CD56 protein
Description	Human recombinant NCAM-1/CD56 protein expressed in Baculovirus with His-Tag. Sequence domain: 20-603aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P13591, 20-603aa; LQVDIVP SQGEISVGES KFFLCQVAGD AKDKDISWFS PNGEKLTNPQ QRISVVWNDD SSSTLTIYNA NIDDAGIYKC VVTGEDGSES EATVNVKIFQ KLMFKNAPT QEFREGEDAV IVCDVVSSLP PTIIWKHKGR DVILKKDVRF IVLSNNYLQI RGIKKTDEGT YRCEGRILAR GEINFKDIQV IVNVPPTIQA RQNIVNATAN LGQSVTLVCD AEGFPEPTMS WTKDGEQIEQ EEDDEKYIFS DDSSQLTIKK VDKNDEAEYI CIAENKAGEQ DATIHLKVFA KPKITYVENQ TAMELEEQVT LTCEASGDPI PSITWRTSTR NISSEKTLT GHMVVRSHAR

VSSLTKSIQ YTDAGEYICT ASNTIGQDSQ SMYLEVQYAP KLQGPVAVYT WEGNQVNITC EVFAYPSATI SWFRDGQLLP SSNYSNIKIY
NTPSASYLEV TPDSENDGYN YNCTAVNRIG QESLEFILVQ ADTPSSPSID QVEPYSSTAQ VQFDEPEATG GVPILKYKAE WRAVGEEVWH
SKWYDAKEAS MEGIVTIVGL KPETTYAVRL AALNGKGLGE ISAASEFKTQ PVHSPPP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human recombinant NCAM-1/CD56 protein

For Research Use Only. Not for use in diagnostic procedures.