

Human CEACAM-1/CD66a Recombinant Protein

Catalog Number: PROTP13688

Overview

Product Name	Human CEACAM-1/CD66a Recombinant Protein
Description	Human CEACAM-1/CD66a Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 35-428aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P13688, 35-428aa; QLTTE SMPFNVAEGK EVLLLVHNL P QQLFGYSWYK GERVDGNRQI VGYAIGTQQA TPGPANS GRE TIYPNASLLI QNVTQNDTGF YTLQVIKSDL VNEEATGQFH VYPELPKPSI SSNNSNPVED KDAVAFTCEP ETQDTTYLWW INNQSLPVSP RLQLSNGNRT LTLLSVTRND TGPYECEIQN PVSANRSDPV TLNVTYGPD TPTISPSDTYY RPGANLSLSC YAASNPPAQY SWLINGTFQQ STQELFIPNI TVNNSGSYTC HANNSVTGCN RTTVKTIIVT ELSPVAKPQ IKASKTTVTG DKDSVNL TCS TNDTGISIRW FFKNQSLPSS ERMKLSQGN T TLSINPVKRE DAGTYWCEVF NPISKNQSDP IMLNVN YNAL PQENGLSPG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CEACAM-1/CD66a Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.