

## Human Enolase 3/Muscle-specific enolase Recombinant Protein

Catalog Number: PROTP13929

### Overview

Product Name	Human Enolase 3/Muscle-specific enolase Recombinant Protein
Description	Human Enolase 3/Muscle-specific enolase Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-434aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% glycerol, 0.1M NaCl

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P13929, 1-434aa; MAMQKIFARE ILDSRGNPTV EVDLHTAKGR FRAAVPSGAS TGIYEALRLR DGDKGRYLGK GVLKAVENIN STLGPALLQK KLSVADQEKV DKFMIELDGT ENKSKFGANA ILGVSLAVCK AGAAEKGVPL YRHIADLAGN PDLILPVPF NVINGGSHAG NKLAMQEFMI LPVGASSFKE AMRIGAEVYH HLKGVKAKY GKDATNVGDE GGFAPNILEN NEALELLKTA IQAAGYPAQ VIGMDVAASE FYRNGKYDLD FKSPDDPARH ITGEKLGELY KSFKNYPVV SIEDPFDQDD WATWTSFLSG VNIQIVGDDL TVTNPKRIAQ AVEKKACNCL LLKVNQIGSV TESIQAACKLA QSNWGWGMVS HRSGETEDTF IADLVVGLCT GQIKTGAPCR SERLAKYNQL MRIEEALGDK AIFAGRKFRN PKAK

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Enolase 3/Muscle-specific enolase Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.