

## Human HSD17B1 Recombinant Protein

Catalog Number: PROTP14061

### Overview

Product Name	Human HSD17B1 Recombinant Protein
Description	Human HSD17B1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-328aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P14061, 1-328aa; MARTVV LITGCSSGIG LHLAVRLASD PSQSFKYYAT LRD LKTQGR L WEAARALACP PGSLETQLD VRDSKSVAAA  
RERVTEGRVD VLVCNAGLGL LGPLEALGED AVASVLDVNV VGTVRMLQAF LPDMKRRGSG RVLVTGVS VGG LMGLPFNDVY CASKFALEGL  
CESLAVLLLP FGVHLSLIEC GPVHTAFMEK VLGSPPEVLD RTDIHTFHRF YQYLAH SKQV FREAAQNPEE VAEVFLTALR APKPTLRYFT  
TERFLPLLRM RLDDPSGSNY VTAMHREVFV DVPKAEAGA EAGGGAGPGA EDEAGR GAVG DPELGDPPAA PQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human HSD17B1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.