

Mouse C-Reactive Recombinant Protein/CRP Recombinant Protein

Catalog Number: PROTP14847

Overview

Product Name	Mouse C-Reactive Recombinant Protein/CRP Recombinant Protein
Description	Mouse C-Reactive Recombinant Protein/CRP Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 20-225aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 0.1M PMSF, 40% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P14847, 20-225aa; HEDMFKKAFV FPKESDTSYV SLEAESKKPL NTFTVCLHFY TALSTVRSFS VFSYATKKNs NDILIFWNKD
KQYTFGVGGA EVRFMVSEIP EAPTHICASW ESATGIVEFW IDGKPKVRKS LHKGYTVGPD ASIILGQEQD SYGGDFDAKQ SLVGDIGDVN
MWDFVLSPEQ ISTVYVGGTL SPNVLNWRAL NYKAQGDVFI KPQLWSLEHH HHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse C-Reactive Recombinant Protein/CRP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.