

Human DARS Recombinant Protein

Catalog Number: PROTP14868

Overview

Product Name	Human DARS Recombinant Protein
Description	Human DARS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-501aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 40% glycerol, 0.15M NaCl, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P14868, 1-501aa; MPSASASRKS QEKPREIMDA AEDYAKERYG ISSMIQSQEK PDRVLVVRVD LTIQKADEVV WVRARVHTSR
AKGKQCFLVL RQQQFNVQAL VAVGDHASKQ MVKFAANINK ESIVDVEGVV RKNVQKIGSC TQQDVELHVQ KIYVISLAEP RLPLQLDDAV
RPEAEGEEEG RATVNQDTRL DNRVIDLRTS TSQAVFRLQS GICHLFRET L INKGFVEIQT PKIISAASEG GANVFTVSYF KNNAYLAQSP
QLYKQMCICA DFEKVFSIGP VFRAEDSNTH RHLTEFVGLD IEMAFNYHYH EVMEEIADTM VQIFKGLQER FQTEIQT VNK QFPCEPFKFL
EPTLRLEYCE ALAMLREAGV EMGDEDDLST PNEKLLGHLV KEKYD TDFYI LDKYPLAVRP FYTMPDPRNP KQSNSYDMFM RGEEILSGAQ
RIHDPQLLTE RALHHGIDLE KIKAYID SFR FGAPPHAGGG IGLERV TMLF LGLHNVRQTS MFPRDPKRLT P

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DARS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.